

## RAC3 抗原（重组蛋白）

中文名称：RAC3 抗原（重组蛋白）

英文名称：RAC3 Antigen (Recombinant Protein)

储 存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 1-192 amino acids of human RAC3

技术规格

<b>Full name:</b>	ras-related C3 botulinum toxin substrate 3 (rho family, small GTP binding protein Rac3)
<b>Swissprot:</b>	P60763
<b>Gene Accession:</b>	BC009605
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Expression system:</b>	Escherichia coli
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Background:</b>	The protein encoded by this gene is a GTPase which belongs to the RAS superfamily of small GTP-binding proteins. Members of this superfamily appear to regulate a diverse array of cellular events, including the control of cell growth, cytoskeletal reorganization, and the activation of protein kinases. Alternative splicing results in multiple transcript variants.