

## 兔抗 IRX4 多克隆抗体

中文名称:	兔抗 IRX4 多克隆抗体
英文名称:	Anti-IRX4 rabbit polyclonal antibody
别名:	iroquois homeobox 4; IRXA3
抗原:	IRX4
储存:	冷冻 (-20℃)
宿主:	Rabbit
反应种属:	Human, Mouse
标记物:	Unconjugate
相关类别:	一抗
克隆类型:	rabbit polyclonal

### 技术规格

<b>Background:</b>	Likely to be an important mediator of ventricular differentiation during cardiac development.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	IRX4
<b>Immunogen:</b>	Synthetic peptide of human IRX4
<b>Full name:</b>	iroquois homeobox 4
<b>Synonyms:</b>	IRXA3
<b>SwissProt:</b>	P78413
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human esophagus cancer and Human lung cancer

<b>IHC Recommend dilution:</b>	40-200
<b>WB Predicted band size:</b>	54 kDa
<b>WB Positive control:</b>	Mouse heart tissue lysate
<b>WB Recommended dilution:</b>	500-2000



