

## 兔抗 IRX5 多克隆抗体

中文名称: 兔抗 IRX5 多克隆抗体

英文名称: Anti-IRX5 rabbit polyclonal antibody

别名: iroquois homeobox 5; HMMS; IRXB2; IRX-2a

抗原: IRX5

储存: 冷冻 (-20°C)

宿主: Rabbit

反应种属: Human, Mouse

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a member of the iroquois homeobox gene family, which are involved in several embryonic developmental processes. Knockout mice lacking this gene show that it is required for retinal cone bipolar cell differentiation, and that it negatively regulates potassium channel gene expression in the heart to ensure coordinated cardiac repolarization. Alternatively spliced transcript variants encoding different isoforms have been found for this gene.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	IRX5
<b>Immunogen:</b>	Synthetic peptide of human IRX5
<b>Full name:</b>	iroquois homeobox 5

<b>Synonyms:</b>	HMMS; IRXB2; IRX-2a
<b>SwissProt:</b>	P78411
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human lung cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	30-150
<b>WB Predicted band size:</b>	50 kDa
<b>WB Positive control:</b>	Mouse lung tissue lysate
<b>WB Recommended dilution:</b>	200-1000





