

## RADIL 抗原（重组蛋白）

中文名称： RADIL 抗原（重组蛋白）

英文名称： RADIL Antigen (Recombinant Protein)

别名： Rap associating with DIL domain; RASIP2

储存： 冷冻（-20℃）

相关类别： 抗原

### 概述

Fusion protein corresponding to a region derived from 876-1075 amino acids of human RADIL

### 技术规格

<b>Background:</b>	Downstream effector of Rap required for cell adhesion and migration of neural crest precursors during development.
<b>Tags:</b>	His tag C-Terminus, GST tag N-Terminus
<b>Expression system:</b>	Escherichia coli
<b>Purity:</b>	>85%, as determined by Coomassie blue stained SDS-PAGE
<b>Gene Accession:</b>	BC004919
<b>Swissprot:</b>	Q96JH8
<b>Synonyms:</b>	RASIP2
<b>Full name:</b>	Rap associating with DIL domain