

兔抗 PARP8 多克隆抗体

中文名称：兔抗 PARP8 多克隆抗体

英文名称：Anti-PARP8 rabbit polyclonal antibody

别名：ARTD16; pART16

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PARP8

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Poly(ADP-ribosylation) is a method of DNA damage-dependent posttranslational modification that helps to rescue injured proliferating cells from cell death. The PARP (poly(ADP-ribose) polymerase) proteins comprise a superfamily of enzymes that functionally modify histones and other nuclear proteins, thereby preventing cell death.
Applications:	ELISA, IHC
Name of antibody:	PARP8
Immunogen:	Synthetic peptide of human PARP8
Full name:	poly (ADP-ribose) polymerase family, member 8
Synonyms :	ARTD16; pART16

SwissProt:	Q8N3A8
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human thyroid cancer and Human brain
IHC Recommend dilution:	50-200

