

兔抗 IRX6 多克隆抗体

中文名称: 兔抗 IRX6 多克隆抗体

英文名称: Anti-IRX6 rabbit polyclonal antibody

别名: iroquois homeobox 6; IRX7; IRX-3; IRXB3

储存: 冷冻 (-20℃)

抗原: IRX6

宿主: Rabbit

反应种属: Human

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	IRX6 is one Iroquois (Irx) proteins comprise a family of homeodomain-containing transcription factors involved in patterning and regionalization of embryonic tissues in both vertebrates and invertebrates.
Applications:	ELISA, IHC
Name of antibody:	IRX6
Immunogen:	Synthetic peptide of human IRX6
Full name:	iroquois homeobox 6
Synonyms:	IRX7; IRX-3; IRXB3
SwissProt:	P78412
ELISA Recommended dilution	5000-10000

n:	
IHC positive control:	Human tonsil
IHC Recommend dilution:	20-100

