

兔抗 SCAF4 多克隆抗体

中文名称：兔抗 SCAF4 多克隆抗体

英文名称： Anti-SCAF4 rabbit polyclonal antibody

别名： SRA4; SFRS15

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SCAF4

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:

This gene likely encodes a member of the arginine/serine-rich splicing factor family. A similar protein in Rat appears to bind the large subunit of RNA polymerase II and provide a link between transcription

	and pre-mRNA splicing. Alternatively spliced transcript variants have been described.
Applications:	WB
Name of antibody:	SCAF4
Immunogen:	Synthesized peptide derived from internal of human SFRS15.
Full name:	SR-related CTD-associated factor 4
Synonyms :	SRA4; SFRS15
SwissProt:	O95104
WB Predicted band size:	126 kDa
WB Positive control:	COLO cells lysate
WB Recommended dilution:	500-3000

