

## 兔抗 PATJ 多克隆抗体

- 中文名称：兔抗 PATJ 多克隆抗体
- 英文名称：Anti-PATJ rabbit polyclonal antibody
- 别名：PATJ, crumbs cell polarity complex component; Cipp; INADL; hINADL; InaD-like
- 相关类别：一抗
- 储存：冷冻（-20℃）
- 宿主：Rabbit
- 抗原：PATJ
- 反应种属：Human
- 标记物：Unconjugate
- 克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a protein with multiple PDZ domains. PDZ domains mediate protein-protein interactions, and proteins with multiple PDZ domains often organize multimeric complexes at the plasma membrane. This protein localizes to tight junctions and to the apical membrane of epithelial cells. A similar protein in Drosophila is a scaffolding protein which tethers several members of a multimeric signaling complex in photoreceptors.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	PATJ
<b>Immunogen:</b>	Synthetic peptide of human PATJ

<b>Full name:</b>	PATJ, crumbs cell polarity complex component
<b>Synonyms:</b>	Cipp; INADL; hINADL; InaD-like
<b>SwissProt:</b>	Q8NI35
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human cervical cancer
<b>IHC Recommend dilution:</b>	40-200

