

RASSF3 抗原（重组蛋白）

中文名称： RASSF3 抗原（重组蛋白）

英文名称： RASSF3 Antigen (Recombinant Protein)

别名： Ras association (RalGDS/AF-6) domain family member 3; RASSF5

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 71-238 amino acids of human RASSF3

技术规格

Full name:	Ras association (RalGDS/AF-6) domain family member 3
Synonyms:	RASSF5
Swissprot:	Q86WH2
Gene Accession:	BC100952
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The RAS oncogene (MIM 190020) is mutated in nearly one-third of all human cancers. Members of the RAS superfamily are plasma membrane GTP-binding proteins that modulate intracellular signal transduction pathways. A subfamily of RAS effectors, including RASSF3, share a RAS association (RA) domain.