

RASSF4 抗原（重组蛋白）

中文名称： RASSF4 抗原（重组蛋白）

英文名称： RASSF4 Antigen (Recombinant Protein)

别名： Ras association (RalGDS/AF-6) domain family member 4; AD037

储存： 冷冻（-20℃）

相关类别： 抗原

概述

Fusion protein corresponding to a region derived from 127-321 amino acids of human RASSF4

技术规格

Full name:	Ras association (RalGDS/AF-6) domain family member 4
Synonyms:	AD037
Swissprot:	Q9H2L5
Gene Accession:	BC032593
Purity:	>85%, as determined by Coomassie blue stained SDS-PAGE
Expression system:	Escherichia coli
Tags:	His tag C-Terminus, GST tag N-Terminus
Background:	The function of this gene has not yet been determined but may involve a role in tumor suppression. Alternative splicing of this gene results in several transcript variants; however, most of the variants have not been fully described.