

## 兔抗 ISM1 多克隆抗体

中文名称：兔抗 ISM1 多克隆抗体

英文名称：Anti-ISM1 rabbit polyclonal antibody

别名：isthmin 1; ISM; Isthmin; C20orf82; bA149I18.1; dJ1077I2.1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ISM1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Acts as an angiogenesis inhibitor.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	ISM1
<b>Immunogen:</b>	Synthetic peptide of human ISM1
<b>Full name:</b>	isthmin 1
<b>Synonyms:</b>	ISM; Isthmin; C20orf82; bA149I18.1; dJ1077I2.1
<b>SwissProt:</b>	B1AKI9
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil
<b>IHC Recommend dilution:</b>	50-300
<b>WB Predicted band size:</b>	52 kDa

<b>WB Positive control:</b>	HepG2 cell lysate
<b>WB Recommended dilution:</b>	1000-5000



