

## 兔抗 SCCPDH 多克隆抗体

中文名称:	兔抗 SCCPDH 多克隆抗体
英文名称:	Anti-SCCPDH rabbit polyclonal antibody
别名:	saccharopine dehydrogenase (putative); NET11; CGI-49
相关类别:	一抗
储存:	冷冻 (-20℃)
宿主:	Rabbit
抗原:	SCCPDH
反应种属:	Human, Mouse, Rat
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

### 技术规格

<b>Background:</b>	SCCPDH belongs to the saccharopine dehydrogenase family. The exact function of SCCPDH remains unknown.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	SCCPDH
<b>Immunogen:</b>	Fusion protein of human SCCPDH
<b>Full name:</b>	saccharopine dehydrogenase (putative)
<b>Synonyms:</b>	NET11; CGI-49
<b>SwissProt:</b>	Q8NBX0
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and Human esophagus cancer
<b>IHC Recommend dilution:</b>	50-300

<b>WB Predicted band size:</b>	47 kDa
<b>WB Positive control:</b>	Human fetal liver tissue, PC-3 and A172 cell lysates
<b>WB Recommended dilution:</b>	1000-5000



