

## 兔抗 ISM2 多克隆抗体

中文名称： 兔抗 ISM2 多克隆抗体

英文名称： Anti-ISM2 rabbit polyclonal antibody

别名： TAIL1; THSD3

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ISM2

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene contains a type 1 thrombospondin domain, which is present in thrombospondin, a number of proteins involved in the complement pathway, as well as in extracellular matrix proteins. Two alternatively spliced transcript variants encoding distinct isoforms have been observed.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	ISM2
<b>Immunogen:</b>	Synthetic peptide of humanISM2
<b>Full name:</b>	isthmin 2
<b>Synonyms :</b>	TAIL1; THSD3

<b>SwissProt:</b>	Q6H9L7
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human liver cancer and Human esophagus cancer
<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	64 kDa
<b>WB Positive control:</b>	293T cells
<b>WB Recommended dilution:</b>	200-1000



