

兔抗 SCG3 多克隆抗体

中文名称：兔抗 SCG3 多克隆抗体

英文名称： Anti-SCG3 rabbit polyclonal antibody

别名： SGIII

储存： 冷冻（-20℃）

抗原： SCG3

宿主： Rabbit

相关类别： 一抗

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a member of the chromogranin/secretogranin family of neuroendocrine secretory proteins. Granins may serve as precursors for biologically active peptides. Some granins have been shown to function as helper proteins in sorting and proteolytic processing of prohormones; however, the function of this protein is unknown. Two transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, WB, IHC
Name of antibody:	SCG3
Immunogen:	Fusion protein of human SCG3

Full name:	Secretogranin III
Synonyms :	SGIII
SwissProt:	Q8WXD2
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human brain
IHC Recommend dilution:	15-50
WB Predicted band size:	53 kDa
WB Positive control:	Mouse intestinum tenue tissue
WB Recommended dilution:	200-1000



