

兔抗 SCGB1A1 多克隆抗体

中文名称：兔抗 SCGB1A1 多克隆抗体

英文名称： Anti-SCGB1A1 rabbit polyclonal antibody

别名： UGB; UP1; CC10; CC16; CCSP; UP-1; CCPBP

储存： 冷冻（-20℃）

抗原： SCGB1A1

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a member of the secretoglobin family of small secreted proteins. The encoded protein has been implicated in numerous functions including anti-inflammation, inhibition of phospholipase A2 and the sequestering of hydrophobic ligands. Defects in this gene are associated with a susceptibility to asthma.
Applications:	ELISA, IHC
Name of antibody:	SCGB1A1
Immunogen:	Fusion protein of human SCGB1A1
Full name:	secretoglobin, family 1A, member 1 (uteroglobin)

Synonyms :	UGB; UP1; CC10; CC16; CCSP; UP-1; CCPBP
SwissProt:	P11684
ELISA Recommended dilution:	2000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	50-200

