

兔抗 PAX9 多克隆抗体

中文名称：兔抗 PAX9 多克隆抗体

英文名称： Anti-PAX9 rabbit polyclonal antibody

别名： STHAG3

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： PAX9

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene is a member of the paired box (PAX) family of transcription factors. Members of this gene family typically contain a paired box domain, an octapeptide, and a paired-type homeodomain. These genes play critical roles during fetal development and cancer growth. The specific function of the paired box 9 gene is unknown but it may involve development of stratified squamous epithelia as well as various organs and skeletal elements.
Applications:	ELISA, IHC
Name of antibody:	PAX9
Immunogen:	Synthetic peptide of human PAX9

Full name:	paired box 9
Synonyms :	STHAG3
SwissProt:	P55771
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human thyroid cancer and Human gastric cancer
IHC Recommend dilution:	25-100

