

兔抗 ISOC2 多克隆抗体

中文名称：兔抗 ISOC2 多克隆抗体

英文名称：Anti-ISOC2 rabbit polyclonal antibody

别名：isochorismatase domain containing 2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ISOC2

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The exact function of ISOC2 remains unknown.
Applications:	ELISA, IHC
Name of antibody:	ISOC2
Immunogen:	Fusion protein of human ISOC2
Full name:	isochorismatase domain containing 2
SwissProt:	Q96AB3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human thyroid cancer
IHC Recommend dilution:	100-300



