

兔抗 PAXIP1 多克隆抗体

中文名称：兔抗 PAXIP1 多克隆抗体

英文名称：Anti-PAXIP1 rabbit polyclonal antibody

别名：PAX interacting protein 1; PTIP; TNRC2; CAGF28; CAGF29; PACIP1; PAXIP1L

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：PAXIP1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene is a member of the paired box (PAX) gene family and encodes a nuclear protein with six BRCT (breast cancer carboxy-terminal) domains. This protein plays a critical role in maintaining genome stability, condensation of chromatin and progression through mitosis.
Applications:	ELISA, IHC
Name of antibody:	PAXIP1
Immunogen:	Fusion protein of human PAXIP1
Full name:	PAX interacting protein 1
Synonyms:	PTIP; TNRC2; CAGF28; CAGF29; PACIP1; PAXIP1L
SwissProt:	Q6ZW49

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human gastric cancer
IHC Recommend dilution:	50-200

