

兔抗 VENTX 多克隆抗体

中文名称：兔抗 VENTX 多克隆抗体

英文名称：Anti-VENTX rabbit polyclonal antibody

别名：NA88A; HPX42B; VENTX2

储存：冷冻（-20℃）

抗原：VENTX

宿主：Rabbit

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

相关类别：一抗

包装：

25 ul 400

100 ul 1500

200 ul 2500

技术规格

Background:	This gene encodes a member of the Vent family of homeodomain proteins. The encoded protein may function as a transcriptional repressor and be involved in mesodermal patterning and hemopoietic stem cell maintenance. Multiple pseudo genes exist for this gene. A transcribed pseudogene located on chromosome X may lead to antigen production in certain melanomas.
Applications:	ELISA, IHC
Name of antibody:	VENTX
Immunogen:	Full length fusion protein

Full name:	VENT homeobox
Synonyms:	NA88A; HPX42B; VENTX2
SwissProt:	O95231
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommended dilution:	20-100