

## 兔抗 SCGB2A2 多克隆抗体

中文名称：兔抗 SCGB2A2 多克隆抗体

英文名称： Anti-SCGB2A2 rabbit polyclonal antibody

别名： MGB1; UGB2

相关类别： 一抗

储存： 冷冻（-20℃） 避光

抗原： SCGB2A2

宿主： Rabbit

反应种属： Human Rat

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

|                          |  |
|--------------------------|--|
| <b>Background:</b>       | Belongs to the secretoglobin family. Lipophilin subfamily. 2 isoforms of the human protein are produced by alternative splicing. |
| <b>Applications:</b>     | WB, IHC, IF  |
| <b>Name of antibody:</b> | SCGB2A2  |

|                                 |   |
|---------------------------------|---|
| <b>Immunogen:</b>               | Synthesized peptide derived from human mam maglobin . |
| <b>Full name:</b>               | secretoglobin, family 2A, member 2                    |
| <b>Synonyms :</b>               | MGB1; UGB2  |
| <b>SwissProt:</b>               | Q13296  |
| <b>IHC positive control:</b>    | Human breast carcinoma tissue                         |
| <b>IHC Recommend dilution:</b>  | 50-100  |
| <b>WB Predicted band size:</b>  | 10 kDa  |
| <b>WB Positive control:</b>     | HepG2 cells lysate                                    |
| <b>WB Recommended dilution:</b> | 500-3000  |
| <b>IF positive control:</b>     | HepG2 cells   |
| <b>IF Recommend dilution:</b>   | 100-500   |



