

兔抗 VEZT 多克隆抗体

中文名称：兔抗 VEZT 多克隆抗体

英文名称：Anti-VEZT rabbit polyclonal antibody

级 别：VEZATIN

抗 原：VEZT

储 存：冷冻（-20℃）

宿 主：Rabbit

相关类别：一抗

反应种属：Human

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a transmembrane protein which has been localized to adherens junctions and shown to bind to myosin VIIA. Examination of expression of this gene in gastric cancer tissues have shown that expression is decreased which appears to be related to hypermethylation of the promoter. Expression of this gene may also be inhibited by binding of a specific microRNA to a target sequence in the 3' UTR of the transcripts. A pseudo gene of this gene is located on the X chromosome. Alternative splicing results in multiple transcript variants.
Applications:	ELISA, IHC
Name of antibody:	VEZT
Immunogen:	Synthetic peptide of human VEZT

Full name:	vezatin, adherens junctions transmembrane protein
Synonyms:	VEZATIN
SwissProt:	Q9HBM0
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer and human tonsil
IHC Recommend dilution:	20-100



