

兔抗 VGLL1 多克隆抗体

中文名称：兔抗 VGLL1 多克隆抗体

英文名称：Anti-VGLL1 rabbit polyclonal antibody

别名：TDU; VGL1

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：VGLL1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene binds proteins of the TEA domain family of transcription factors (TEFs) through the Vg (vestigial) homology region found in its N-terminus. It may thus function as a specific coactivator for the mammalian TEFs.
Applications:	ELISA, WB
Name of antibody:	VGLL1
Immunogen:	Synthetic peptide of human VGLL1
Full name:	vestigial like family member 1
Synonyms:	TDU; VGL1

SwissProt:	Q99990
ELISA Recommended dilution:	2000-20000
WB Predicted band size:	29 kDa
WB Positive control:	293T cells
WB Recommended dilution:	2000-10000

