

## 兔抗 SCGB3A1 多克隆抗体

中文名称：兔抗 SCGB3A1 多克隆抗体

英文名称：Anti-SCGB3A1 rabbit polyclonal antibody

级 别：HIN1; HIN-1; LU105; UGRP2; PnSP-2

抗 原：SCGB3A1

储 存：冷冻（-20℃）

宿 主：Rabbit

相关类别：一抗

反应种属：Human, Mouse

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Potential growth inhibitory cytokine.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SCGB3A1
<b>Immunogen:</b>	Synthetic peptide of human SCGB3A1
<b>Full name:</b>	secretoglobin family 3A member 1
<b>Synonyms:</b>	HIN1; HIN-1; LU105; UGRP2; PnSP-2
<b>SwissProt:</b>	Q96QR1
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human liver cancer
<b>IHC Recommend dilution:</b>	25-100

