

## 兔抗 SCGN 多克隆抗体

中文名称：兔抗 SCGN 多克隆抗体

英文名称： Anti-SCGN rabbit polyclonal antibody

别名： secretagoin, EF-hand calcium binding protein; SEGN; CALBL; SECRET; setagin; DJ501N12.8

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SCGN

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	The encoded protein is a secreted calcium-binding protein which is found in the cytoplasm. It is related to calbindin D-28K and calretinin. This protein is thought to be involved in KCL-stimulated calcium flux and cell proliferation.
<b>Applications:</b>	ELISA, WB
<b>Name of antibody:</b>	SCGN
<b>Immunogen:</b>	Fusion protein of human SCGN
<b>Full name:</b>	secretagoin, EF-hand calcium binding protein
<b>Synonyms:</b>	SEGN; CALBL; SECRET; setagin; DJ501N12.8

<b>SwissProt:</b>	O76038
<b>ELISA Recommended dilution:</b>	5000-10000
<b>WB Predicted band size:</b>	32 kDa
<b>WB Positive control:</b>	Hela cell and Human fetal liver tissue lysates
<b>WB Recommended dilution:</b>	500-2000

