

兔抗 VGLL2 多克隆抗体

中文名称：兔抗 VGLL2 多克隆抗体

英文名称：Anti-VGLL2 rabbit polyclonal antibody

别名：VGL2; VITO1

相关类别：一抗

储存：冷冻（-20℃）

抗原：VGLL2

宿主：Rabbit

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	This gene encodes a protein with a transcriptional enhancer factor 1 (TEF-1) interaction domain. The encoded protein may act as a co-factor of TEF-1 regulated gene expression during skeletal muscle development. Alternatively spliced transcript variants encoding multiple isoforms have been observed.
Applications:	ELISA, IHC
Name of antibody:	VGLL2
Immunogen:	Full length fusion protein
Full name:	vestigial like family member 2
Synonyms:	VGL2; VITO1

SwissProt:	Q8N8G2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human colorectal cancer
IHC Recommend dilution:	25-100



