

兔抗 SCMH1 多克隆抗体

中文名称：兔抗 SCMH1 多克隆抗体

英文名称：Anti-SCMH1 rabbit polyclonal antibody

别名：Scm polycomb group protein homolog 1; Scml3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SCMH1

反应种属：Human, Mouse

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	Associates with Polycomb group (PcG) multiprotein complexes; the complex class is required to maintain the transcriptionally repressive state of some genes.
Applications:	ELISA, WB, IHC
Name of antibody:	SCMH1
Immunogen:	Fusion protein of human SCMH1
Full name:	Scm polycomb group protein homolog 1
Synonyms:	Scml3
SwissProt:	Q96GD3
ELISA Recommended dilution:	5000-10000

IHC positive control:	Human cervical cancer
IHC Recommend dilution:	100-300
WB Predicted band size:	73 kDa
WB Positive control:	Raji and 231 cell lysates
WB Recommended dilution:	500-2000



