

兔抗 VHL 多克隆抗体

中文名称：兔抗 VHL 多克隆抗体

英文名称： Anti-VHL rabbit polyclonal antibody

别名： RCA1; VHL1; pVHL; HRCA1

抗原： VHL

储存： 冷冻（-20℃）

宿主： Rabbit

相关类别： 一抗

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Von Hippel-Lindau syndrome (VHL) is a dominantly inherited familial cancer syndrome predisposing to a variety of malignant and benign tumors. A germline mutation of this gene is the basis of familial inheritance of VHL syndrome. The protein encoded by this gene is a component of the protein complex that includes elongin B, elongin C, and cullin-2, and possesses ubiquitin ligase E3 activity.
Applications:	ELISA, WB, IHC
Name of antibody:	VHL
Immunogen:	Full length fusion protein

Full name:	von Hippel-Lindau tumor suppressor, E3 ubiquitin protein ligase
Synonyms :	RCA1; VHL1; pVHL; HRCA1
SwissProt:	P40337
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human colon cancer and Human thyroid cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	24 kDa
WB Positive control:	Human fetal brain tissue
WB Recommended dilution:	500-2000



