

兔抗 SCN2B 多克隆抗体

中文名称：兔抗 SCN2B 多克隆抗体

英文名称： Anti-SCN2B rabbit polyclonal antibody

别名： ATFB14

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SCN2B

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:

Crucial in the assembly, expression, and functional modulation of the heterotrimeric complex of the sodium channel. The subunit beta-2 causes an increase in the plasma membrane surface area and in its folding into mic

	rovilli. Interacts with TNR may play a crucial role in clustering and regulation of activity of sodium channels at nodes of Ranvier By similarity.
Applications:	WB
Name of antibody:	SCN2B
Immunogen:	Synthesized peptide derived from internal of human SCN2B.
Full name:	sodium channel, voltage-gated, type II, beta subunit
Synonyms :	ATFB14
SwissProt:	O60939
WB Predicted band size:	24 kDa
WB Positive control:	293 cells lysate
WB Recommended dilution :	500-3000

