

兔抗 VIL1 多克隆抗体

中文名称: 兔抗 VIL1 多克隆抗体

英文名称: Anti-VIL1 rabbit polyclonal antibody

别名: VIL; D2S1471

抗原: VIL1

储存: 冷冻 (-20°C)

宿主: Rabbit

相关类别: 一抗

反应种属: Human

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a member of a family of calcium-regulated actin-binding proteins. This protein represents a dominant part of the brush border cytoskeleton which functions in the capping, severing, and bundling of actin filaments. Two mRNAs of 2.7 kb and 3.5 kb have been observed; they result from utilization of alternate poly-adenylation signals present in the terminal exon.
Applications:	ELISA, WB, IHC
Name of antibody:	VIL1
Immunogen:	Full length fusion protein
Full name:	villin 1

Synonyms :	VIL; D2S1471
SwissProt:	P09327
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human prostate cancer and Human esophagus cancer
IHC Recommend dilution:	25-100
WB Predicted band size:	93 kDa
WB Positive control:	HT-29 cell
WB Recommended dilution:	200-1000



