

## 兔抗 PCK2 多克隆抗体

中文名称：兔抗 PCK2 多克隆抗体

英文名称： Anti-PCK2 rabbit polyclonal antibody

别名： PEPCK; PEPCK2; PEPCK-M

相关类别： 一抗

抗原： PCK2

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human, Mouse

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	This gene encodes a mitochondrial enzyme that catalyzes the conversion of oxaloacetate to phosphoenolpyruvate in the presence of guanosine triphosphate (GTP). A cytosolic form of this protein is encoded by a different gene and is the key enzyme of gluconeogenesis in the liver. Alternatively spliced transcript variants have been described.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PCK2
<b>Immunogen:</b>	Fusion protein of human PCK2
<b>Full name:</b>	phosphoenolpyruvate carboxykinase 2
<b>Synonyms :</b>	PEPCK; PEPCK2; PEPCK-M

<b>SwissProt:</b>	Q16822
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human brain and human liver cancer
<b>IHC Recommend dilution:</b>	50-200
<b>WB Predicted band size:</b>	71 kDa
<b>WB Positive control:</b>	293T, Jurkat and Hela cells
<b>WB Recommended dilution:</b>	500-2000



