

兔抗 SCN4B 多克隆抗体

中文名称：兔抗 SCN4B 多克隆抗体

英文名称： Anti-SCN4B rabbit polyclonal antibody

别名： Lag; SMN; OP18; PP17; PP19; PR22; LAP18; C1orf215

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： SCN4B

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:

This gene belongs to the stathmin family of genes. It encodes a ubiquitous cytosolic phosphoprotein proposed to function as an intracellular relay integrating regulatory signals of the cellular e

	<p>nvironment. The encoded protein is involved in the regulation of the microtubule filament system by destabilizing microtubules. It prevents assembly and promotes disassembly of microtubules. Multiple transcript variants encoding different isoforms have been found for this gene.</p>
Applications:	WB
Name of antibody:	SCN4B
Immunogen:	Synthesized peptide derived from internal of human SCN4B.
Full name:	stathmin 1
Synonyms :	Lag; SMN; OP18; PP17; PP19; PR22; LA P18; C1orf215
SwissProt:	P16949
WB Predicted band size:	17 kDa
WB Positive control:	HepG2 cells and Jurkat cells lysates
WB Recommended dilution:	500-3000

