

## 兔抗 VIM(Phospho-Ser56) 多克隆抗体

中文名称：兔抗 VIM(Phospho-Ser56) 多克隆抗体

英文名称： Anti-VIM(Phospho-Ser56) rabbit polyclonal antibody

别名： HEL113; CTRCT30

相关类别： 一抗

储存： 冷冻（-20℃） 避光

宿主： Rabbit

抗原： VIM(Phospho-Ser56)

反应种属： Human

标记物： Unconjugate

克隆类型： Unconjugate

### 技术规格

<b>Background:</b>	This gene encodes a member of the intermediate filament family. Intermediate filaments, along with microtubules and actin microfilaments, make up the cytoskeleton. The protein encoded by this gene is responsible for maintaining cell shape, integrity of the cytoplasm, and stabilizing cytoskeletal interactions. It is also involved in the immune response, and controls the transport of low-density lipoprotein (LDL)-derived cholesterol from a lysosome to the site of esterification.
<b>Applications:</b>	WB

<b>Name of antibody:</b>	VIM(Phospho-Ser56)
<b>Immunogen:</b>	Peptide sequence around phosphorylation site of serine 56 (A-S-S(p)-P-G) derived from Human Vimentin.
<b>Full name:</b>	vimentin
<b>Synonyms :</b>	HEL113; CTRCT30
<b>SwissProt:</b>	P08670
<b>WB Predicted band size:</b>	54 kDa
<b>WB Positive control:</b>	A549 cells lysates treated with Nocodazole
<b>WB Recommended dilution:</b>	500-1000

