

兔抗 VIPR1 多克隆抗体

中文名称：兔抗 VIPR1 多克隆抗体

英文名称： Anti-VIPR1 rabbit polyclonal antibody

别名： II; HVR1; RDC1; V1RG; VIPR; VIRG; VAPC1; VPAC1; VPAC1R; VIP-R-1; VPCAP1R; PACAP-R2; PACAP-R-2

相关类别： 一抗

储存： 冷冻 (-20℃)

宿主： Rabbit

抗原： VIPR1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes a receptor for vasoactive intestinal peptide, a small neuropeptide. Vasoactive intestinal peptide is involved in smooth muscle relaxation, exocrine and endocrine secretion, and water and ion flux in lung and intestinal epithelia. Its actions are effected through integral membrane receptors associated with a guanine nucleotide binding protein which activates adenylate cyclase. Several transcript variants encoding different isoforms have been found for this gene.
Applications:	ELISA, IHC

Name of antibody:	VIPR1
Immunogen:	Synthetic peptide of humanVIPR1
Full name:	vasoactive intestinal peptide receptor 1
Synonyms :	II; HVR1; RDC1; V1RG; VIPR; VIRG; VAPC1; VPAC1; VPA C1R; VIP-R-1; VPCAP1R; PACAP-R2; PACAP-R-2
SwissProt:	P32241
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human liver cancer and Human brain
IHC Recommend dilution:	10-50

