

兔抗 ITGA7 多克隆抗体

中文名称：兔抗 ITGA7 多克隆抗体

英文名称：Anti-ITGA7 rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：ITGA7

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	<p>The protein encoded by this gene belongs to the integrin alpha chain family. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. They mediate a wide spectrum of cell-cell and cell-matrix interactions, and thus play a role in cell migration, morphological development, differentiation, and metastasis. This protein functions as a receptor for the basement membrane protein laminin-1. It is mainly expressed in skeletal and cardiac muscles and may be involved in differentiation and migration processes during myogenesis. Defects in this gene are associated with congenital myopathy. Alternatively spliced transcript variants encoding different isoforms have been noted for this gene.</p>
--------------------	---

Applications:	ELISA, IHC
Name of antibody:	ITGA7
Immunogen:	Synthetic peptide of human ITGA7
Full name:	integrin subunit alpha 7
SwissProt:	Q13683
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human brain
IHC Recommend dilution:	25-100

