

兔抗 VIP 多克隆抗体

中文名称: 兔抗 VIP 多克隆抗体

英文名称: Anti-VIP rabbit polyclonal antibody

别名: PHM27

抗原: VIP

储存: 冷冻 (-20℃)

宿主: Rabbit

反应种属: Human, Rat

相关类别: 一抗

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	The protein encoded by this gene belongs to the glucagon family. It stimulates myocardial contractility, causes vaso dilation, increases glycogenolysis, lowers arterial blood pressure and relaxes the smooth muscle of trachea, stomach and gall bladder. The protein also acts as an antimicrobial peptide with antibacterial and antifungal activity. Alternative splicing occurs at this locus and two transcript variants encoding distinct isoforms have been identified.
Applications:	ELISA, IHC
Name of antibody:	VIP
Immunogen:	Full length fusion protein

Full name:	vasoactive intestinal peptide
Synonyms :	PHM27
SwissProt:	P01282
ELISA Recommended dilution:	2000-5000
IHC positive control:	Human liver cancer and Human prostate cancer
IHC Recommend dilution:	25-100

