

## 兔抗 SCN10A 多克隆抗体

中文名称：兔抗 SCN10A 多克隆抗体

英文名称：Anti-SCN10A rabbit polyclonal antibody

别名：PN3; SNS; hPN3; Nav1.8

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SCN10A

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	The protein encoded by this gene is a tetrodotoxin-resistant voltage-gated sodium channel alpha subunit. The properties of the channel formed by the encoded transmembrane protein can be altered by interaction with different beta subunits. This protein may be involved in the onset of pain associated with peripheral neuropathy.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SCN10A
<b>Immunogen:</b>	Synthetic peptide of human SCN10A
<b>Full name:</b>	sodium channel, voltage-gated, type X, alpha subunit
<b>Synonyms :</b>	PN3; SNS; hPN3; Nav1.8

<b>SwissProt:</b>	Q9Y5Y9
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human breast cancer and Human cervical cancer
<b>IHC Recommend dilution:</b>	15-50

