

兔抗 VNN1 多克隆抗体

中文名称：兔抗 VNN1 多克隆抗体

英文名称：Anti-VNN1 rabbit polyclonal antibody

别名：HDLCQ8; Tiff66

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：VNN1

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

This gene encodes a member of the vanin family of proteins, which share extensive sequence similarity with each other, and also with biotinidase. The family includes secreted and membrane-associated proteins, a few of which have been reported to participate in hematopoietic cell trafficking. No biotinidase activity has been demonstrated for any of the vanin proteins, however, they possess pantetheinase activity, which may play a role in oxidative-stress response. This protein, like its mouse homolog, is likely a GPI-anchored cell surface molecule. The mouse protein is expressed by the perivascular thymic stroma

	I cells and regulates migration of T-cell progenitors to the thymus. This gene lies in close proximity to, and in the same transcriptional orientation as, two other vanin genes on chromosome 6q23-q24.
Applications:	ELISA, IF
Name of antibody:	VNN1
Immunogen:	Synthetic peptide of human VNN1
Full name:	vanin 1
Synonyms :	HDLCQ8; Tiff66
SwissProt:	O95497
ELISA Recommended dilution:	5000-240000
IF positive control:	HepG2 and NCCIT cell
IF Recommended dilution:	100-500

