

## 兔抗 PCTP 多克隆抗体

中文名称:	兔抗 PCTP 多克隆抗体
英文名称:	Anti-PCTP rabbit polyclonal antibody
别名:	phosphatidylcholine transfer protein; PC-TP; STARD2
抗原:	PCTP
储存:	冷冻 (-20℃)
宿主:	Rabbit
反应种属:	Human, Mouse
相关类别:	一抗
标记物:	Unconjugate
克隆类型:	rabbit polyclonal

### 技术规格

<b>Background:</b>	Catalyzes the transfer of phosphatidylcholine between membranes. Binds a single lipid molecule.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	PCTP
<b>Immunogen:</b>	Fusion protein of human PCTP
<b>Full name:</b>	phosphatidylcholine transfer protein
<b>Synonyms:</b>	PC-TP; STARD2
<b>SwissProt:</b>	Q9UKL6
<b>ELISA Recommended dilution:</b>	5000-10000
<b>IHC positive control:</b>	Human tonsil and Human esophagus cancer
<b>IHC Recommend dilution:</b>	40-200

<b>WB Predicted band size:</b>	25 kDa
<b>WB Positive control:</b>	Human fetal liver tissue, HepG2 cell, Human placenta tissue lysates
<b>WB Recommended dilution:</b>	200-1000



