

兔抗 ITGAL 多克隆抗体

中文名称：兔抗 ITGAL 多克隆抗体

英文名称：Anti-ITGAL rabbit polyclonal antibody

别名：CD11A, LFA-1, LFA1A

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：ITGAL

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

ITGAL encodes the integrin alpha L chain. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. This I-domain containing alpha integrin combines with the beta 2 chain (ITGB2) to form the integrin lymphocyte function-associated antigen-1 (LFA-1), which is expressed on all leukocytes. LFA-1 plays a central role in leukocyte intercellular adhesion through interactions with its ligands, ICAMs 1-3 (intercellular adhesion molecules 1 through 3), and also functions in lymphocyte costimulatory signaling. Two transcript variants encoding different isoforms have been found for this gene.

Applications:	ELISA, IHC
Name of antibody:	ITGAL
Immunogen:	Synthetic peptide of human ITGAL
Full name:	integrin, alpha L (antigen CD11A (p180), lymphocyte function-associated antigen 1; alpha polypeptide)
Synonyms :	CD11A, LFA-1, LFA1A
SwissProt:	P20701
ELISA Recommended dilution:	1000-2000
IHC positive control:	Human esophagus cancer and Human liver cancer
IHC Recommend dilution:	15-50

