

兔抗 VPREB3 多克隆抗体

中文名称：兔抗 VPREB3 多克隆抗体

英文名称：Anti-VPREB3 rabbit polyclonal antibody

别名：8HS20; N27C7-2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：VPREB3

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

The protein encoded by this gene is the human ortholog of the mouse VpreB3 (8HS20) protein, is thought to be involved in B-cell maturation, and may play a role in assembly of the pre-B cell receptor (pre-BCR). While the role of this protein in B-cell development has not yet been elucidated, studies with the chicken ortholog of this protein have found that when overexpressed, this protein localizes to the endoplasmic reticulum. The mouse ortholog of this protein has been shown to associate with membrane mu heavy chains early in the course of pre-B cell receptor biosynthesis

	. Expression of this gene has been observed in some lymphomas.
Applications:	ELISA, IHC
Name of antibody:	VPREB3
Immunogen:	Fusion protein of human VPRES3
Full name:	pre-B lymphocyte 3
Synonyms:	8HS20; N27C7-2
SwissProt:	Q9UKI3
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human cervical cancer
IHC Recommend dilution:	25-100



