

兔抗 ITGAX 多克隆抗体

中文名称： 兔抗 ITGAX 多克隆抗体

英文名称： Anti-ITGAX rabbit polyclonal antibody

别名： CD11C, SLEB6

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： ITGAX

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	This gene encodes the integrin alpha X chain protein. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. This protein combines with the beta 2 chain (ITGB2) to form a leukocyte-specific integrin referred to as inactivated-C3b (iC3b) receptor 4 (CR4). The alpha X beta 2 complex seems to overlap the properties of the alpha M beta 2 integrin in the adherence of neutrophils and monocytes to stimulated endothelium cells, and in the phagocytosis of complement coated particles.
Applications:	ELISA, IHC
Name of antibody:	ITGAX

Immunogen:	Synthetic peptide of human ITGAX
Full name:	integrin, alpha X (complement component 3 receptor 4 subunit)
Synonyms:	CD11C, SLEB6
SwissProt:	P20702
IHC positive control:	Human gastric cancer
IHC Recommend dilution :	15-50

