

兔抗 PDCD1LG2 多克隆抗体

中文名称：兔抗 PDCD1LG2 多克隆抗体

英文名称： Anti-PDCD1LG2 rabbit polyclonal antibody

别名： B7DC; Btdc; PDL2; CD273; PD-L2; PDCD1L2; bA574F11.2

相关类别： 一抗

抗原： PDCD1LG2

储存： 冷冻（-20℃）

宿主： Rabbit

反应种属： Human

标记物： Unconjugate

克隆类型： rabbit polyclonal

技术规格

Background:	Involved in the costimulatory signal, essential for T-cell proliferation and IFNG production in a PDCD1-independent manner. Interaction with PDCD1 inhibits T-cell proliferation by blocking cell cycle progression and cytokine production (By similarity).
Applications:	ELISA, IHC
Name of antibody:	PDCD1LG2
Immunogen:	Fusion protein of human PDCD1LG2
Full name:	programmed cell death 1 ligand 2
Synonyms:	B7DC; Btdc; PDL2; CD273; PD-L2; PDCD1L2; bA574F11.2
SwissProt:	Q9BQ51

ELISA Recommended dilution:	5000-10000
IHC positive control:	Human liver cancer and Human esophagus cancer
IHC Recommend dilution:	50-200



