

## 兔抗 SCPEP1 多克隆抗体

中文名称：兔抗 SCPEP1 多克隆抗体

英文名称： Anti-SCPEP1 rabbit polyclonal antibody

别名： RISC; HSCP1

相关类别： 一抗

储存： 冷冻（-20℃）

宿主： Rabbit

抗原： SCPEP1

反应种属： Human, Mouse, Rat

标记物： Unconjugate

克隆类型： rabbit polyclonal

### 技术规格

<b>Background:</b>	Retinoid-inducible serine carboxypeptidase is an enzyme that in humans is encoded by the SCPEP1 gene. May be involved in vascular wall and kidney homeostasis.
<b>Applications:</b>	ELISA, WB, IHC
<b>Name of antibody:</b>	SCPEP1
<b>Immunogen:</b>	Fusion protein of human SCPEP1
<b>Full name:</b>	serine carboxypeptidase 1
<b>Synonyms :</b>	RISC; HSCP1
<b>SwissProt:</b>	Q9HB40
<b>ELISA Recommended dilution:</b>	1000-2000
<b>IHC positive control:</b>	Human cervical cancer

<b>IHC Recommend dilution:</b>	25-100
<b>WB Predicted band size:</b>	51 kDa
<b>WB Positive control:</b>	231 cells
<b>WB Recommended dilution:</b>	200-1000



