

兔抗 SCRN2 多克隆抗体

中文名称：兔抗 SCRN2 多克隆抗体

英文名称：Anti-SCRN2 rabbit polyclonal antibody

别名：secernin 2; Ses2

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SCRN2

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:

SCRN genes encode protein members of the secernin family including Secernin 1, Secernin 2 and Secernin 3. Function of Secernin 2 and Secernin 3 are so far not well understood. Secernin 1 is a novel 50-kDa cytosolic protein that appears to be involved in the regulation of exocytosis from peritoneal mast cells. Secernin-1 a novel tumor-associated antigen (TAA) and may be a universal marker of different cancer types including gastric cancer, which has been validated by several researches. Measurement of SCRNI can be used in the early detecti

	on of cancer or in the surveillance of patients who undergo surgery.
Applications:	ELISA, WB, IHC
Name of antibody:	SCRN2
Immunogen:	Fusion protein of human SCRIN2
Full name:	secernin 2
Synonyms:	Ses2
SwissProt:	Q96FV2
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human esophagus cancer and Human tonsil
IHC Recommend dilution:	50-300
WB Predicted band size:	47 kDa
WB Positive control:	Mouse small intestines tissue lysate
WB Recommended dilution:	500-2000



