

兔抗 SCTR 多克隆抗体

中文名称：兔抗 SCTR 多克隆抗体

英文名称：Anti-SCTR rabbit polyclonal antibody

别名：SR

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SCTR

反应种属：Human

标记物：Unconjugate

克隆类型：rabbit polyclonal

技术规格

Background:	The protein encoded by this gene is a G protein-coupled receptor and belongs to the glucagon-VIP-secretin receptor family. It binds secretin which is the most potent regulator of pancreatic bicarbonate, electrolyte and volume secretion. Secretin and its receptor are suggested to be involved in pancreatic cancer and autism.
Applications:	ELISA, WB, IHC
Name of antibody:	SCTR
Immunogen:	Synthetic peptide of human SCTR
Full name:	secretin receptor

Synonyms:	SR
SwissProt:	P47872
ELISA Recommended dilution:	5000-10000
IHC positive control:	Human tonsil
IHC Recommend dilution:	20-100
WB Predicted band size:	50 kDa
WB Positive control:	Human stomach tissue lysate
WB Recommended dilution:	200-1000



