

兔抗 VPS25 多克隆抗体

中文名称: 兔抗 VPS25 多克隆抗体

英文名称: Anti-VPS25 rabbit polyclonal antibody

别名: DERP9; EAP20; FAP20

抗原: VPS25

储存: 冷冻 (-20℃)

宿主: Rabbit

相关类别: 一抗

反应种属: Human, Mouse, Rat

标记物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a protein that is a subunit of the endosomal sorting complex required for transport II (ESCRT-II). This protein complex functions in sorting of ubiquitinated membrane proteins during endocytosis. A pseudogene of this gene is present on chromosome 1.
Applications:	ELISA, WB
Name of antibody:	VPS25
Immunogen:	Full length fusion protein
Full name:	vacuolar protein sorting 25 homolog

Synonyms:	DERP9; EAP20; FAP20
SwissProt:	Q9BRG1
ELISA Recommended dilution:	2000-5000
WB Predicted band size:	21 kDa
WB Positive control:	Human fetal liver tissue
WB Recommended dilution:	500-2000

