

## 兔抗 SCT 多克隆抗体

中文名称：兔抗 SCT 多克隆抗体

英文名称：Anti-SCT rabbit polyclonal antibody

相关类别：一抗

储 存：冷冻（-20℃）

宿 主：Rabbit

抗 原：SCT

反应种属：Human, Mouse, Rat

标 记 物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

<b>Background:</b>	Secretin belongs to the glucagon family. This protein is an endocrine hormone and its major site of production is the endocrine S cells located in the proximal small intestinal mucosa. The release of active secretin is stimulated by either fatty acids or an acidic pH in the duodenum. This hormone stimulates the secretion of bicarbonate-rich pancreatic fluids and has also been shown to regulate the growth and development of the stomach, small intestine, and pancreas. Secretin deficiency has been implicated in autistic syndrome, suggesting that the hormone could have a neuroendocrine function in addition to its role in digestion.
<b>Applications:</b>	ELISA, IHC

<b>Name of antibody:</b>	SCT
<b>Immunogen:</b>	Synthetic peptide of human SCT
<b>Full name:</b>	secretin
<b>SwissProt:</b>	P09683
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human thyroid cancer and Human ovarian cancer
<b>IHC Recommend dilution:</b>	50-200

