

兔抗 PDCD4(Phospho-Ser67) 多克隆抗体

中文名称：兔抗 PDCD4(Phospho-Ser67) 多克隆抗体

英文名称： Anti-PDCD4(Phospho-Ser67) rabbit polyclonal antibody

别名： H731

抗原： PDCD4(Phospho-Ser67)

储存： 冷冻（-20℃） 避光

宿主： Rabbit

相关类别： 一抗

反应种属： Human Mouse Rat

标记物： Unconjugate

克隆类型： Unconjugate

技术规格

Background:	This gene is a tumor suppressor and encodes a protein that binds to the eukaryotic translation initiation factor 4A1 and inhibits its function by preventing RNA binding. Alternative splicing results in multiple transcript variants.
Applications:	WB
Name of antibody:	PDCD4(Phospho-Ser67)
Immunogen:	Peptide sequence around phosphorylation site of serine 67 (K-N-S(p)-S-R) derived from Human PDCD4 .
Full name:	programmed cell death 4 (neoplastic transformation inhibitor)
Synonyms :	H731

SwissProt:	Q53EL6
WB Predicted band size:	52 kDa
WB Positive control:	HepG2 cells lysates
WB Recommended dilution:	500-1000

