

兔抗 VPS28 多克隆抗体

中文名称: 兔抗 VPS28 多克隆抗体

英文名称: Anti-VPS28 rabbit polyclonal antibody

储 存: 冷冻 (-20℃)

抗 原: VPS28

宿 主: Rabbit

反应种属: Human, Mouse, Rat

相关类别: 一抗

标 记 物: Unconjugate

克隆类型: rabbit polyclonal

技术规格

Background:	This gene encodes a protein subunit of the ESCRT-I complex (endosomal complexes required for transport), which functions in the transport and sorting of proteins into subcellular vesicles. This complex can also be hijacked to facilitate the budding of enveloped viruses from the cell membrane. Alternative splicing results in multiple transcript variants encoding different isoforms.
Applications:	ELISA, IHC
Name of antibody:	VPS28
Immunogen:	Full length fusion protein
Full name:	VPS28 subunit of ESCRT-I
SwissProt:	Q9UK41
ELISA Recommended dilution:	2000-5000

IHC positive control:	Human thyroid cancer and Human prostate cancer
IHC Recommend dilution:	25-100



