

## 兔抗 SDC3 多克隆抗体

中文名称：兔抗 SDC3 多克隆抗体

英文名称：Anti-SDC3 rabbit polyclonal antibody

别名：SDCN; SYND3

相关类别：一抗

储存：冷冻（-20℃）

宿主：Rabbit

抗原：SDC3

反应种属：Human, Mouse, Rat

标记物：Unconjugate

克隆类型：rabbit polyclonal

### 技术规格

#### Background:

The protein encoded by this gene belongs to the syndecan proteoglycan family. It may play a role in the organization of cell shape by affecting the actin cytoskeleton, possibly by transferring signals from the cell surface

	in a sugar-dependent mechanism. Allelic variants of this gene have been associated with obesity.
<b>Applications:</b>	ELISA, IHC
<b>Name of antibody:</b>	SDC3
<b>Immunogen:</b>	Synthetic peptide of human SDC3
<b>Full name:</b>	syndecan 3
<b>Synonyms :</b>	SDCN; SYND3
<b>SwissProt:</b>	O75056
<b>ELISA Recommended dilution:</b>	2000-5000
<b>IHC positive control:</b>	Human brain
<b>IHC Recommend dilution:</b>	40-150

